



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Electric Service Rubber Gloves

**Product Code :**  
LBNY-0002-660008



**Description :**

Electric Service Rubber Gloves

**Technical Specification :**

Made for 11" Rubber Insulated Gloves.  
Manufactured from top grain goat skin.  
Provide excellent protection for rubber insulating gloves.  
All Salisbury Leather protectors  
Length: 10".

Electric Service Rubber Gloves, Electric Service Rubber Gloves Equipments, Electric Service Rubber Gloves Tools, Electric Service Rubber Gloves Tool Kits, Electric Service Rubber Gloves Manufacturers, Electric Service Rubber Gloves Suppliers from India, China, Kenya



**Laboratory Instrument India**