



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Laboratory Steam Generator

**Product Code :**  
LBNY-0005-1100007



**Description :**

Easy and safe access for maintenance and inspection, with door isolator to the main electrical enclosure.  
A floor-standing continuous steam generator.  
Solid state level control on feed water tank.  
Splash proof electrical enclosure and switches.  
Manufactured in accordance with Module G of the European Pressure Equipment directive.  
A comprehensive instruction manual.

**Technical Specification :**

Laboratory Steam Generator



**Laboratory Instrument India**