



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laboratory Steam Generator

Product Code :
LBNY-0005-1100007



Description :

Easy and safe access for maintenance and inspection, with door isolator to the main electrical enclosure.
A floor-standing continuous steam generator.
Solid state level control on feed water tank.
Splash proof electrical enclosure and switches.
Manufactured in accordance with Module G of the European Pressure Equipment directive.
A comprehensive instruction manual.

Technical Specification :

Laboratory Steam Generator



Laboratory Instrument India