



Email : sales@laboratoryinstrumentindia.com

Product Name :
Corrosion Studies Kit

Product Code :
LBNY-0005-11000016



Description :

The corrosion studies kit allows simultaneous study of up to eight corrosion cells each containing three test specimens.

Corrosion rates measured visually and by weighing.

Stirring by air or inert gas.

Sample mounting in a manner to minimise secondary effects.

Comprising:

A service panel

Air pump

Test cell supports

Platinum electrodes

Low voltage supply

Digital pH meter

Specimen pieces of steel, zinc, brass and copper

Associated glassware, tubing and wiring.

Eight test cells of glass construction, with specially machined lids enabling samples to be mounted.

Technical Specification :

Corrosion Studies Kit



Laboratory Instrument India