



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Leather Protectors Glove

**Product Code :**  
LBNY-0002-660009



**Description :**

Leather Protectors Glove

**Technical Specification :**

Leather protector gloves always be worn over Insulating Rubber Gloves to provide mechanical protection against cuts, abrasions, and punctures.

Made from top grain goat leather, the protectors increase the life of your rubber gloves.

Length: 10".

Leather Protectors Glove, Leather Protectors Glove Equipments, Leather Protectors Glove Tools, Leather Protectors Glove Tool Kits, Leather Protectors Glove Manufacturers, Leather Protectors Glove Suppliers from India, China, Kenya



**Laboratory Instrument India**