



Email : sales@laboratoryinstrumentindia.com

Product Name :
Leather Protectors Glove

Product Code :
LBNY-0002-660009



Description :

Leather Protectors Glove

Technical Specification :

Leather protector gloves always be worn over Insulating Rubber Gloves to provide mechanical protection against cuts, abrasions, and punctures.

Made from top grain goat leather, the protectors increase the life of your rubber gloves.

Length: 10".

Leather Protectors Glove, Leather Protectors Glove Equipments, Leather Protectors Glove Tools, Leather Protectors Glove Tool Kits, Leather Protectors Glove Manufacturers, Leather Protectors Glove Suppliers from India, China, Kenya



Laboratory Instrument India