



Email : sales@laboratoryinstrumentindia.com

Product Name :
Solid Liquid Extraction Unit

Product Code :
LBNY-0005-1100019



Description :

Self contained floor standing unit with integral control and instrumentation console
Probes at each stage allow full monitoring of the process With the data logger, the product concentrations can be directly displayed in real time (for the recommended material)
Three stage extraction process, can also be configured as one or two stages
Independent vessel for closed loop batch extraction
Temperature control at each stage
Individual control of solvent feed pumps at each stage, of cell rotation rate and of feed rate
Optional data logging accessory with USB interface and educational software.

Technical Specification :

Solid Liquid Extraction Unit



Laboratory Instrument India