



Email : sales@laboratoryinstrumentindia.com

Product Name :
EV Battery Lift (Capacity 1 Ton)

Product Code :
LBNY-0002-660010



Description :

EV Battery Lift (Capacity 1 Ton)

Technical Specification :

Easy to move with 4 high strength universal self-locking casters
Providing the height limit to ensure safe operation
Durable self-lubricating sliding blocks prolong the service life of the machine
The platform can move to 25mm in any direction horizontally
The platform can be tilted in vertical directions
Carefully select premium quality parts and components
Fully covered with insulating rubber to ensure safe and non-slip operation.

EV Battery Lift (Capacity 1 Ton), EV Battery Lift (Capacity 1 Ton) Equipments, EV Battery Lift (Capacity 1 Ton) Tools, EV Battery Lift (Capacity 1 Ton) Tool Kits, EV Battery Lift (Capacity 1 Ton) Manufacturers, EV Battery Lift (Capacity 1 Ton) Suppliers from India, China, Kenya



Laboratory Instrument India

