



Email : sales@laboratoryinstrumentindia.com

Product Name :
Water Chilled Recirculating Unit

Product Code :
LBNY-0005-11200012



Description :

Refrigeration unit to chill water:
Power consumption: 370 W.
Refrigerating capacity: 40W at -10/32°C.
Expansion valve.
Cooling compressor.
Evaporator and condenser.
Coolant accumulation tank.
Tank of division of the cooling liquid.
Cooling filter.

Technical Specification :

Centrifugal pump:
Water tank, volume: 15 L.
Max. head: 30 m.
Max. flow rate: 600 L/h.
Power consumption: 120 W.
Enthalpy diagram of the refrigerant R134a.
Control panel:
Electronic controller for water temperature.
Cables and Accessories, for normal operation.
Unit ON/OFF switch.
Manuals:

Dimensions: 400 x 250 x 500 mm. approx. (15.74 x 9.84 x 19.68 inches approx.)

Weight: 20 Kg.

Water Chilled Recirculating Unit is supplied with the following manuals: Required Services, Assembly and Installation, Starting-up, Safety, Maintenance & Practices Manuals.



Laboratory Instrument India