




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Water Chilled Recirculating Unit

**Product Code :**  
LBNY-0005-11200012



**Description :**

Refrigeration unit to chill water:  
Power consumption: 370 W.  
Refrigerating capacity: 40W at -10/32°C.  
Expansion valve.  
Cooling compressor.  
Evaporator and condenser.  
Coolant accumulation tank.  
Tank of division of the cooling liquid.  
Cooling filter.

**Technical Specification :**

Centrifugal pump:  
Water tank, volume: 15 L.  
Max. head: 30 m.  
Max. flow rate: 600 L/h.  
Power consumption: 120 W.  
Enthalpy diagram of the refrigerant R134a.  
Control panel:  
Electronic controller for water temperature.  
Cables and Accessories, for normal operation.  
Unit ON/OFF switch.  
Manuals:

---

Dimensions: 400 x 250 x 500 mm. approx. (15.74 x 9.84 x 19.68 inches approx.)

Weight: 20 Kg.

Water Chilled Recirculating Unit is supplied with the following manuals: Required Services, Assembly and Installation, Starting-up, Safety, Maintenance & Practices Manuals.



**Laboratory Instrument India**