



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Insulated Rescue Hooks (6ft, 1000V)

**Product Code :**  
LBNY-0002-660011



**Description :**

Insulated Rescue Hooks (6ft, 1000V)

**Technical Specification :**

Insulated rescue hook is an invaluable tool for any workplace.

It is used to release a victim from contact.

Confined spaces, in vaults, or just near electrical cabinets and switch gear are some of the places where this tool is a must.

Insulated rescue hooks feature a foam-filled, fiberglass reinforced handle for superior electrical insulation and a coated, heat treated body hook with an 18" opening and are.

Insulated Rescue Hooks (6ft, 1000V), Insulated Rescue Hooks (6ft, 1000V) Equipments, Insulated Rescue Hooks (6ft, 1000V) Tools, Insulated Rescue Hooks (6ft, 1000V) Tool Kits, Insulated Rescue Hooks (6ft, 1000V) Manufacturers, Insulated Rescue Hooks (6ft, 1000V) Suppliers from India, China, Kenya



**Laboratory Instrument India**

