




Email : sales@laboratoryinstrumentindia.com

Product Name :
Automobile Air Conditioning Trainer

Product Code :
LBNY-0005-11200014



Description :

Automobile Air Conditioning Trainer

Bench-top unit.

Diagram in the front panel with similar distribution to the elements in the real unit.

A condenser (radiator).

Anodized aluminum structure and panels of painted steel.

Main metallic elements in stainless steel.

An evaporator.

Two computer controlled refrigeration fans.

A filter.

A computer controlled variable area compressor.

An expansion valve.

A storage tank.

A computer controlled electric motor (to simulate the cars engine).

A variable frequency drive to adjust the number of revolutions of the motor between 300 and 1500 r.p.m.

A safety pressure switch that stops the compressor at 14 bars.

A compartment made of methacrylate with ventilation windows.

Technical Specification :

Automobile Air Conditioning Trainer



Laboratory Instrument India