



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Computer Controlled Cooling Plant With Ice Store

**Product Code :**  
LBNY-0005-11200015



#### **Description :**

Cooling plant at teaching and industrial level, with modular design, and computer controlled.

Plant with ice store, wet cooling tower and dry cooler.

Main metallic elements in stainless steel.

Diagram in the front panel.

The different units (modules) connected with hoses.

Anodized aluminium structures and panels in painted steel.

Wet cooling tower.

Dry cooler.

Using valves different operating modes can be configured.

Refrigeration circuit (condenser, evaporator, compressor) and pumps.

Manometers.

Temperature sensors.

Flow sensors.

Ice tank.

Liquid tank.

Wattmeter.

Refrigerant R134a.

Connecting hoses.

High pressure control.

Pressure sensors.

Adjustment elements.

#### **Technical Specification :**

---

Computer Controlled Cooling Plant With Ice Store



**Laboratory Instrument India**