



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Refrigeration And Air Conditioning Modular Unit

**Product Code :**  
LBNY-0005-11200019



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Anodized aluminum structure and panels in painted steel.  
Diagram in the front panel with similar distribution to the elements in the real unit  
Main metallic elements in stainless steel.  
Unit to perform experiments with a SCADA SOFTWARE.  
System control via solenoid valves.  
The unit includes wheels to facilitate its mobility.  
Base unit of the unit that allows the development basic experiments on the operation of the refrigeration and air conditioning systems in conjunction with different models.  
Function educational.  
Models attached securely with fasteners.  
Condenser unit with: Compressor, Condenser, Receiver of refrigerant.  
Data acquisition.  
Automation of the air operation system.

**Technical Specification :**

Refrigeration And Air Conditioning Modular Unit

---

# Laboratory Instrument India

