



Email : sales@laboratoryinstrumentindia.com

Product Name :
Heat Pump Air Conditioning Refrigeration Unit With Cycle
Inversion Valve and Two Evaporators

Product Code :
LBNY-0005-11200022



Description :

Bench-top unit, with anodized aluminum structure and painted steel panels.
Main metallic elements of stainless steel.
Water condenser, capacity: 1580 W.
High pressure control.
Coolant accumulator.
Diagram in the front panel with similar distribution to the elements in the real unit.
Cooling compressor, computer controlled, power: 1/2 CV.
Air condenser, computer controlled.
Coolant filter to retain particles of condensate.
Expansion valve.
4 Manometers, situated at important points of the unit.
Water evaporator, capacity: 1580 W.
Air evaporator, computer controlled.
Liquid separator to retain liquid particles before going on to the compressor.

Technical Specification :

Heat Pump Air Conditioning Refrigeration Unit With Cycle Inversion Valve and Two Evaporators

Laboratory Instrument India

