



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Heat Pump Air Conditioning Refrigeration Unit With Cycle  
Inversion Valve and Two Evaporators

**Product Code :**  
LBNY-0005-11200022



**Description :**

Bench-top unit, with anodized aluminum structure and painted steel panels.  
Main metallic elements of stainless steel.  
Water condenser, capacity: 1580 W.  
High pressure control.  
Coolant accumulator.  
Diagram in the front panel with similar distribution to the elements in the real unit.  
Cooling compressor, computer controlled, power: 1/2 CV.  
Air condenser, computer controlled.  
Coolant filter to retain particles of condensate.  
Expansion valve.  
4 Manometers, situated at important points of the unit.  
Water evaporator, capacity: 1580 W.  
Air evaporator, computer controlled.  
Liquid separator to retain liquid particles before going on to the compressor.

**Technical Specification :**

Heat Pump Air Conditioning Refrigeration Unit With Cycle Inversion Valve and Two Evaporators

---

# Laboratory Instrument India

