




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Computer Controlled Contact Plate Freezer

**Product Code :**  
LBNY-0005-11200024



**Description :**

Bench-top unit.  
Main metallic elements in stainless steel.  
Diagram in the front panel with similar distribution to the elements in the real unit.  
Anodized aluminum structure and panels of painted steel.  
Air condenser.  
High pressure control.  
Coolant accumulation tank.  
2 Manometers.  
Coolant compressor, computer controlled.  
Evaporator-freezer, with two freezing plates of 180 mm x 280 mm. One freezing plate is fixed and the other one is height adjustable.  
Four-way valve.  
Expansion valve.  
6 Temperature sensors, type J (temperature measurement of the food).  
8 Temperature sensors:  
2 Temperature sensors, type J (temperature measurement of the coolant).

**Technical Specification :**

Computer Controlled Contact Plate Freezer



**Laboratory Instrument India**