



Email : sales@laboratoryinstrumentindia.com

Product Name :
Instantaneous Electric Boiler for Steam Production

Product Code :
LBNY-0005-11200032



Description :

The water feeding can be directly connected to the mains and the water circulation is ensured by a pump rated at 0.55 kW. This electric and automatic steam generator is fitted with two separate sets of heating elements for single or double operation.

Features:

Steam production: 50 kg/h

Operating pressure: 4-5 bar.

Total electric power: 40 kW at 380 W, three-phase, 50 Hz.

Technical Specification :

Weight and Dimensions:

Dimensions: 830 x 830 x 570 h mm

Weight: 113 kg.



Laboratory Instrument India

