



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Instantaneous Electric Boiler for Steam Production

**Product Code :**  
LBNY-0005-11200034



**Description :**

The water feeding can be directly connected to the mains and the water circulation is ensured by a pump rated at 0.55 kW. This electric and automatic steam generator is fitted with two separate sets of heating elements for single or double operation.

**Technical Specification :**

Steam production: 50 kg/h  
Operating pressure: 4-5 bar.  
Total electric power: 40 kW at 380 W, three-phase, 50 Hz.

Weight and dimensions:  
Dimensions: 830 x 830 x 570 h mm  
Weight: 113 kg.



**Laboratory Instrument India**