



Email : sales@laboratoryinstrumentindia.com

Product Name :
Instantaneous Electric Boiler for Steam Production

Product Code :
LBNY-0005-11200034



Description :

The water feeding can be directly connected to the mains and the water circulation is ensured by a pump rated at 0.55 kW. This electric and automatic steam generator is fitted with two separate sets of heating elements for single or double operation.

Technical Specification :

Steam production: 50 kg/h
Operating pressure: 4-5 bar.
Total electric power: 40 kW at 380 W, three-phase, 50 Hz.

Weight and dimensions:
Dimensions: 830 x 830 x 570 h mm
Weight: 113 kg.



Laboratory Instrument India