



Email : sales@laboratoryinstrumentindia.com

Product Name :
Exhaust Gas Analyzer For Internal Combustion Engines

Product Code :
LBNY-0005-11200039



Description :

Composition and Technical Characteristics:

Combustion Engines

Basic unit complete with:

Battery charger housing.

LAN Ethernet connector.

Infrared ray input by remote control.

Keyboard, mouse, Bluetooth via USB flash drive, serial and USB ports.

Rev counter and universal thermometer with data output both via cable and via Bluetooth, complete with connecting cables to acquire data from the battery and the EOBD diagnostic socket.

Engine oil temperature probe: 5 to 180°C.

PC module, monitor, A4 printer and DVD/CD reader.

User manual, software.

Weight and Dimensions:

Overall dimensions: 760 x 660 x 1570 mm

Weight: 86 kg.

The Exhaust Gas Analyzer For Internal Combustion Engines measures the carbon monoxide (CO), carbon dioxide (CO₂), unburnt hydrocarbons (HC) and oxygen (O₂) values and uses them to determine the air-fuel ratio correction coefficient automatically it is also used to measure smoke emission from diesel engines.

Technical Specification :

Exhaust Gas Analyzer For Internal Combustion Engines



Laboratory Instrument India