



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Computerized Heat Exchanger Bench Top Study Unit

**Product Code :**  
LBNY-0005-11200042



**Description :**

A bench top service unit makes it possible to supply hot water to different types of heat exchangers which are small scale versions of the commonly used industrial ones. The Computerized Heat Exchanger Bench Top Study Unit has been designed to carry out experiments concerning the main techniques of heat transfer. Thus the student can compare their performance, fully understand the importance of a correct choice and, if possible, compare the influence of co-current and counter-current heat exchange.

**Technical Specification :**

Computerized Heat Exchanger Bench Top Study Unit



**Laboratory Instrument India**