



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Computerized Automotive Engine Test Bed

**Product Code :**  
LBNY-0005-11200047



**Description :**

The system includes a dynamometer with a control and instrumentation panel. On the test bed, an engine selected among the engines in widespread use is mounted. Computerized Automotive Engine Test Bed allows to carry out a vast selection of tests on petrol and diesel engines, with a power output of up to 80 kW. Other engines not listed in this catalogue can be supplied on request. The system is supplied with manuals, which describe all parts of the trainer, the installation and utilization procedures, as well as many exercises with the relative results.

**Technical Specification :**

Required Services:  
Water discharge.  
Exhaust gas output.  
Electrical power supply: 220 V single phase, 50/60 Hz, 1 kW;  
Water supply: 5000l/h, 3 bar.



**Laboratory Instrument India**

