



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pensky Martens Flash Point Apparatus

**Product Code :**  
LBNY-0005-11200048



**Description :**

The Pensky Martens Flash Point Apparatus allows for the determination of the flash point up to 370 °C of drying oils, fuel oils, lube oils, suspensions of solids, solvent-type waxes, liquids that tend to form a surface film and other liquids.

**Technical Specification :**

**Requirements Services:**

Gas or LPG supply

Electrical supply: 220V A.C. single-phase, 50/60 Hz, 700 VA

Weight and Dimensions:

Dimensions: 300 x 300 x 600 h mm.

Weight: 10 kg.



**Laboratory Instrument India**