



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pensky Martens Flash Point Apparatus

Product Code :
LBNY-0005-11200048



Description :

The Pensky Martens Flash Point Apparatus allows for the determination of the flash point up to 370 °C of drying oils, fuel oils, lube oils, suspensions of solids, solvent-type waxes, liquids that tend to form a surface film and other liquids.

Technical Specification :

Requirements Services:

Gas or LPG supply

Electrical supply: 220V A.C. single-phase, 50/60 Hz, 700 VA

Weight and Dimensions:

Dimensions: 300 x 300 x 600 h mm.

Weight: 10 kg.



Laboratory Instrument India