



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Junkers Boys Calorimeter

**Product Code :**  
LBNY-0005-11200055



**Description :**

The calorific value is obtained by utilizing a heat exchanger where the products of combustion are cooled by a measured flow of water. Junkers Boys Calorimeter apparatus allows for the accurate determination of the calorific value of the main fuel gases which are used in the industry, such as butane, propane and methane. The system is supplied with manuals which describe all parts of the trainer, the installation and utilization procedures, as well as some exercises with the relative results. The burner is supplied with a complete set of interchangeable nozzles for the different types of gases.

**Technical Specification :**

Weight and Dimensions:  
Dimensions: 1000 x 500 x 1000 h mm  
Weight: 25 kg.



**Laboratory Instrument India**