



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Versatile Data Acquisition System

**Product Code :**  
LBNY-0005-11200065



### Description :

#### Features:

- Display of real-time data, in digital form or as an analogue meter
- Real-time traces of analogue signals
- Logging data for printing and later analysis
- Works with existing Versatile Data Acquisition System compatible products
- Performing real-time calculations to generate user defined data
- Creating and printing charts and data tables
- Recording data manually or automatically
- Data capture set by time or intervals
- Exporting data for use by other software.

#### Technical Specification :

##### Bench Top Interface Unit:

Nett dimensions and weight: Width 330 mm x depth 205 mm x height 45 mm and 2 kg

Packed volume and weight: 0.05 m<sup>3</sup> and 9 kg.

##### Frame-Mounting Interface Unit:

Nett dimensions and weight: Width 190 mm x depth 130 mm x height 445 mm and 4.5 kg

Packed volume and weight: 0.07 m<sup>3</sup> and 9 kg.

---

# Laboratory Instrument India

