



Email : sales@laboratoryinstrumentindia.com

Product Name :
Split Type Air Conditioner System

Product Code :
LBNY-0005-11200067



Description :

Individual components are clearly shown. Temperatures and pressures at various points are indicated. Split Type Air Conditioner System is used for the demonstration of split type room air conditioning operation.

Experiments:

Observation of pressures and temperatures
Study of layout and process diagram
Remote control operation.

Technical Specification :

Technical Data:

Compressor type : Rotary.
Refrigerant : R22.
Cooling capacity : 2.5kW. (approx 9000 BTU).
Measuring instruments
Temperatures : 5 ea.
Remote control : 1 ea.
Pressure gauges : 2 ea.
Power supply : 220V 1Ph 50Hz.

Laboratory Instrument India

