



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Split Type Air Conditioner System

**Product Code :**  
LBNY-0005-11200067



**Description :**

Individual components are clearly shown. Temperatures and pressures at various points are indicated. Split Type Air Conditioner System is used for the demonstration of split type room air conditioning operation.

**Experiments:**

Observation of pressures and temperatures  
Study of layout and process diagram  
Remote control operation.

**Technical Specification :**

**Technical Data:**

Compressor type : Rotary.  
Refrigerant : R22.  
Cooling capacity : 2.5kW. (approx 9000 BTU).  
Measuring instruments  
Temperatures : 5 ea.  
Remote control : 1 ea.  
Pressure gauges : 2 ea.  
Power supply : 220V 1Ph 50Hz.

---

# Laboratory Instrument India

