



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Ice Maker Test Rig

**Product Code :**  
LBNY-0005-11200068



**Description :**

The unit comprises a typical refrigeration system whereby a small ice maker is used for evaporator. Ice Maker Test Rig table top refrigeration system is intended for studying ice making process. Instruments are provided for measurement of pressures and temperatures.

**Experiments:**

Coefficient of performance.  
Familiarization with ice making system and components.

**Technical Specification :**

**Technical Data:**

Ice making unit : 0.4 kg/hr.  
Measuring instruments :  
Temperature digital display.  
Pressure gauges.  
Power supply : 220V 1Ph 50Hz.



**Laboratory Instrument India**

