



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Fiber-Optic Transmission Training System

**Product Code :**  
LBNY-0005-1130004



**Description :**

**Features:**

The equipment that you assemble will transmit voice from one point to another, using light traveling through an optical fiber.

The experiment of shows you how easy it is to make productive use of fiber optic materials.

With four different data transmission ways (self module

Transmission, and module-to-PC transmission).

Transmission, module-to-module transmission, PC-to-module.

**Technical Specification :**

Fiber-Optic Transmission Training System



**Laboratory Instrument India**