



Email : sales@laboratoryinstrumentindia.com

Product Name :
Floor Standing Belt Conveyor

Product Code :
LBNY-0005-11400015



Description :

The Floor Standing Belt Conveyor Learning System adds to the Mechanical Drives 1 Learning System to teach skills related to flat belt conveyors. This system uses industrial quality components to help students become better prepared for what they will encounter on the job and to withstand frequent student use. The model Floor Standing Belt Conveyor includes an operational floor-standing flat belt conveyor; student learning materials for both theory and lab, and teachers guide. Students will learn industry-relevant skills related to this new topic including operation, installation, alignment, maintenance, and troubleshooting. Flat belt conveyors are used in almost every industry to transport material and packages at low cost from one point to another.

Features:

Industry-relevant skills
Teaches skills related to flat belt conveyors.

Technical Specification :

Floor Standing Belt Conveyor



Laboratory Instrument India

