



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Rigging Learning System

**Product Code :**  
LBNY-0005-11400017



**Description :**

The Rigging 2 Learning System expands on the initial offering of rigging knowledge to include synthetic slings, equipment movement, and industrial cranes. Rigging 2 Learning System adds to Rigging 1 to offer a comprehensive learning experience to help learners prepare for certifications. Proper rigging techniques are vital for efficient movement of loads and workers safety.

The Rigging 2 Learning System includes two synthetic slings, four load protection pads, jack, lever dolly, and much more. In addition to equipment and printed curriculum, the offers a student reference guide, instructors manual, installation guide, and an optional online eBook. These components will be used to explore topics such as synthetic sling operation, equipment movement using roller bars and dollies, and crane operation signals.

**Features:**

Replication of Common Loads  
Industrial Standard Components.

**Technical Specification :**

Rigging Learning System



**Laboratory Instrument India**

