



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Centrifugal Pump With Stuffing Box Learning System

**Product Code :**  
LBNY-0005-11400019



#### **Description :**

The Centrifugal Pump With Stuffing Box Learning System includes a cast iron-housed centrifugal pump with stuffing box and a piping network that features a pressure line, suction line, and pressure gauge hose. The Centrifugal Pump With Stuffing Box Learning System explains the installation, function, and maintenance of a centrifugal pump with a stuffing box. The Centrifugal Pump With Stuffing Box Learning System covers a pump design that, in spite of the modern mechanical seal, remains widely used in industry due to its cost-effectiveness and ability to avoid catastrophic failure. Additional topics include: how to troubleshoot this pump design; what the flow/pressure characteristics of a centrifugal pump are; and how to startup and operate a centrifugal pump with a stuffing box. This feature makes it an ideal design for fire departments in fire pump applications. This commitment to providing learners with the most realistic and applicable training available is one example of how is transforming the global workforce one life at a time. These components allow learners to practice skills on an industrial-grade centrifugal pump that they may work with at real-world job sites.

#### **Features:**

Industrial Quality Centrifugal Pump  
Stuffing Box Packing Type Seal  
Mounting Compatible with Workstation.

#### **Technical Specification :**

Centrifugal Pump With Stuffing Box Learning System



**Laboratory Instrument India**