




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Programmable Logic Controller Trainer

**Product Code :**  
LBNY-0005-1150003



#### **Description :**

Input-simulation switches function as level and pulse input for different input signal.  
With various peripheral devices and other devices that support external extensions, particularly suited for laboratory experiment and project implementation  
Installation of output relay helps to increase load current.  
Easy-to-use windows-based development software  
Easy to carry, move and store with a suitcase design  
Source mode (standard) or sink mode (option) upon request.  
Equipped with various simulations I/O devices for the convenience of studying and observing the results  
Use 4mm safety sockets input/output terminals to ensure the physical safety of users.

#### **Technical Specification :**

Programmable Logic Controller Trainer



**Laboratory Instrument India**