



Email : sales@laboratoryinstrumentindia.com

Product Name :
Programmable Logic Controller Trainer PLC

Product Code :
LBNY-0005-1150004



Description :

This trainer contains advanced networking arrangement, distributed control section and process controlling system. It can efficiently regulate hydraulic or pneumatic cylinder. Programmable Logic Controller Trainer PLC has extensive application in automobile, refrigeration, steel manufacturing, textile machine production and glass industries. It is offered as per customers preference on expansion module and PLC type. PC based operation has contributed to its user friendly design. Designed by skilled technocrats, this trainer PLC serves as the reliable option to comprehend complex theories of programmable logic controller. On board IR switches is one of its main factors. We are a prestigious manufacturer and exporter of excellent grade Programmable Logic Controller Trainer PLC. We are a distinguished supplier and manufacturer of Programmable Logic Controller Trainer PLC.

Features:

- Advanced design
- Convenient to maintain
- Easy to understand mechanism
- Simple to handle.

Technical Specification :

Programmable Logic Controller Trainer PLC

Laboratory Instrument India

