



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mechatronics Training System

Product Code :
LBNY-0005-1150005



www.laboratoryinstrumentindia.com

Description :

The system is universal, industry-based, modular and flexible for further expansion. Students can learn the entire process of production such as feeding, processing, etc. Each station simplifies the training of operation and can be expanded sequentially step by step through building complex automated procedure. The modular production system stations allow varying simulation of real production processes that exist in industry field. With 6-Stations modular production system, mechatronics training includes distributing, testing, processing, handling, assembly and storing stations that can be operated separately or integrally.

Technical Specification :

Mechatronics Training System



Laboratory Instrument India