




Email : sales@laboratoryinstrumentindia.com

Product Name :
Sedimentation Tank

Product Code :
LBNY-0005-1160003



Description :

The Sedimentation Tank allows to study the hydraulic and efficiency characteristics of the sedimentation by means of a special rectangular section tank. The unit is fitted with an ink injecting device in order to visualize the flow lines obtained in the settling tank. The sludge, which is kept in suspension in the preparation tank by a stirrer, is sent to the sedimentation tank by a variable speed peristaltic pump together with a water flow rate adjusted through a manual valve. The group is supplied with manuals which describe all parts of the unit, the installation and utilization procedures, as well as many exercises with the relative results.

Technical Specification :

Water supply from mains: 300 l/h
Electrical power supply: 220/240 V single phase, 50/60 Hz; 0.2 kW
Weight and dimensions:
Dimensions: 1500 x 600 x 1400 h mm
Weight: 70 kg.



Laboratory Instrument India

