



Email : sales@laboratoryinstrumentindia.com

Product Name :
Batch Fluidized Bed Drier

Product Code :
LBNY-0005-1160002



Description :

The Batch Fluidized Bed Drier is a small but technologically advanced fluidized bed drier. It allows to directly obtain the characteristic drying curve (CDC), which is used to preset precise parameters on industrial drying plants also of different type such as: cabinets, conveyor plants, rotary. It is suitable to carry out research and small batch production in processes on fluidized materials such as drying, granulation, coatings, etc.

Technical Specification :

Technical Features:

Evaporation capacity: 1.5 2 kg/h H₂O

Feeding capacity: 6 kg

Power: 2.8 kW

Electrical power supply: 220 V CA, 50/60 Hz.

Weight and Dimensions:

Dimensions: 900 x 500 x 1100h mm

Weight: 60 kg.



Laboratory Instrument India

