



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Computerized Multipurpose Distillation Plant

**Product Code :**  
LBNY-0005-1160005



**Description :**

The Computerized Multipurpose Distillation Plant can be used to perform a wide range of tests on binary, continuous or batch distillation processes. A data acquisition and analysis system for Personal Computer, makes it possible to perform the measurements in automatic mode, to produce charts of the characteristics curves which can be displayed on screen or printed, to save the experimental data on the disk or print them out. Suitable options make it possible to process azeotropic mixtures, to carry out the distillation under vacuum or in steam current (stripping). Of modular design, the training set is very flexible and easy to use, and can be integrated at a later stage with one or more of the many options available. A manual of exercises completes the supply.

**Technical Specification :**

Water drainage.

Water supply: 200 l/h - 1 bar minimum

Electrical power supply 220/240 V A.C. single-phase - 50/60 Hz - 3.5 kW.



**Laboratory Instrument India**

