



Email : sales@laboratoryinstrumentindia.com

Product Name :
Microscope-Inverted

Product Code :
LBNY-0001-510001



Description :

Microscope-Inverted

Technical Specification :

Eyepiece: 10X magnification

Objectives: 10X,20X, 40X magnification

Binocular head

Operating Environment: Operating Temperature: +15°C to +30°C

Utility Requirement: Electrical power supply: 220VAC \pm 10%, Resettable over-current mains fuse to be incorporated

Microscope-Inverted, Microscope-Inverted Equipments, Microscope-Inverted Tools, Microscope-Inverted Tool Kits, Microscope-Inverted Manufacturers, Microscope-Inverted Suppliers from India, China, Kenya



Laboratory Instrument India