



Email : sales@laboratoryinstrumentindia.com

Product Name :
Press Filter

Product Code :
LBNY-0005-1160008



Description :

This type of plant is particularly suitable to filter both milk and fruit juices and it can be generally employed for any liquid-solid mixture (not soluble).

Exercises of Press Filter

- Filtering efficiency as a function of treated product turbidity
- Introduction in a pilot line for producing pasteurized milk or juices
- Filtering efficiency as a function of delivery pressure
- Filtering efficiency as a function of feeding flow rate

Technical Specification :

- Electrical supply: 220/380 V three-phase, 50/60 Hz; 0.6 kW
- Weight and dimensions of Press Filter
- Plant overall dimensions: 2000 x 2000 x 1000 h mm
- Total plant weight: 160 kg



Laboratory Instrument India

