



Email : sales@laboratoryinstrumentindia.com

Product Name :
Multipurpose Distillation Plant

Product Code :
LBNY-0005-1160009



Description :

The Multipurpose Distillation Plant can be used to perform a wide range of tests on binary, continuous or batch distillation processes. A manual of exercises completes the supply. Suitable options make it possible to process a zeotropic mixtures, to carry out the distillation under vacuum or in steam current (stripping). Of modular design, the training set is very flexible and easy to use, and can be integrated at a later stage with one or more of the many options available.

Technical Specification :

Water supply: 200 l/h - 1 bar

Electrical power supply 220/240 V A.C. single-phase - 50/60 Hz - 3.5 kW

Water drainage

Weights and Dimensions (basic system) of Multipurpose Distillation Plant

Dimensions: 1600 x 630 x 2900h m

Weight: 260 kg.



Laboratory Instrument India

