



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tubular Flow Reactor

Product Code :
LBNY-0005-11600010



Description :

The chemical tubular reactor is widely used in the industrial practice because it is easy to be made and has good conversion yield that can be got by recirculating the mixture. The unit has been designed to realise an absolutely turbulent flow inside the reactor: two reagents are continuously fed and the obtained reaction is homogenous and approximately isothermal. The Tubular Flow Reactor teaching unit allows to study the reactions which can occur in a tubular reactor using the tubular flow method.

Technical Specification :

Electrical supply: 220/240 V, single phase, 50/60Hz
Weight and dimensions of Tubular Flow Reactor
Dimensions: 1000 x 800 x 1000 h mm
Weight: 70 kg.



Laboratory Instrument India