



Email : sales@laboratoryinstrumentindia.com

Product Name :
Double Effect Concentrator

Product Code :
LBNY-0005-11600013



Description :

The Double Effect Concentrator Plant entirely made of stainless material in order to guarantee a long duration, can be utilized for concentration both at atmospheric pressure and under vacuum. Double Effect Concentrator Plant allows to carry out a complete study on the concentration and to produce small quantities of concentrate for laboratory use.

Technical Specification :

Weight and Dimensions:
Dimensions (mm) 2600x1700x 2100h
Weight (kg)400.



Laboratory Instrument India