



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Double Effect Concentrator

**Product Code :**  
LBNY-0005-11600013



**Description :**

The Double Effect Concentrator Plant entirely made of stainless material in order to guarantee a long duration, can be utilized for concentration both at atmospheric pressure and under vacuum. Double Effect Concentrator Plant allows to carry out a complete study on the concentration and to produce small quantities of concentrate for laboratory use.

**Technical Specification :**

Weight and Dimensions:  
Dimensions (mm) 2600x1700x 2100h  
Weight (kg)400.



**Laboratory Instrument India**