



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Top Bench Modular Reactor

**Product Code :**  
LBNY-0005-11600016



#### **Description :**

The Top Bench Modular Reactor which allows to carry out the exercises on different type of reactors. Then the student can compare the performance of the different reactors and fully understand the importance of choosing the most suitable reactor according to the needed result. The Basic Unit provides water at constant temperature to perform the exercises and allows dosing the reagents within the available reactors.

The Basic Unit can be used along with the following optional reactors:

Batch Reactor  
Cascade Stirred Reactors  
Stirred Tank Continuous Reactor  
Tubular Reactor

#### **Technical Specification :**

Electrical supply: 220V single phase, 50/60 Hz, 1.5 kW  
Weights and dimensions:  
Dimensions: 800 x 400 x 500h mm  
Weight: 50 kg.

---

## Laboratory Instrument India

