



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Wetted Wall Absorption Pilot Plant

**Product Code :**  
LBNY-0005-11600021



**Description :**

Feed tank with a capacity of 45 litres.  
Deoxygenation column manufactured from Plexiglas with a height of 1,600 mm  
Pump for absorption column.  
Pump for deoxygenator.  
Wetted wall column of height 900 mm and inside diameter 34 mm.  
Air compressor.  
2x industrial oxymeters 0 to 20 mg per litre or 0 to 200% with a resolution of mg/l or 0.1% SAT.  
2x industrial measuring sensors for dissolved oxygen.  
Variable area flowmeter for water with glass tube and a range of 1.9 to 19 litres per hour.  
Variable area flowmeter for air with glass tube and a range of 30 to 300 litres per hour.  
Emergency stop button  
Electrical control panel in main cabinet.  
Power Supply: 230 Vac 50 Hz single-phase 1500 VA.

**Technical Specification :**

Wetted Wall Absorption Pilot Plant

---

# Laboratory Instrument India

