



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cathodic Corrosion Protection Training Plant

Product Code :
LBNY-0005-11600023



Description :

Safety glasses and gloves
Digital voltmeter
PC interface for the measurement and data storage
Digital ammeter on console
Digital voltmeter on console
Zn reference cell
Copper electrode
Carbon steel electrode
Cu/CuSO₄ reference cell
Ag/AgCl reference cell
Transparent basin to build the electrolytic test bath
Simple circuit with sliding variable resistor and lamp suitable for inserting into the electrical circuit of the electrolytic cell
Zinc electrode
Aluminum electrode
Magnesium electrode
DC feeder (with constant voltage, constant current, constant potential facilities)
Ti-Pt anode
MMO tubular anode
Fe bar
Cu bar
Fe-Ni bar
Waterproof resistor with thermostatic device

Resistivity fluid cell
Air pump with relevant sprayer
Carbon steel electrode (partially coated with epoxy compound)
Carbon steel electrode (completely coated with epoxy compound)
Various reagents supplied in plastic cans with safety technical data sheets
Set of ancillaries and connecting leads
Set of spare fuses.

Technical Specification :

Cathodic Corrosion Protection Training Plant



Laboratory Instrument India