



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cathodic Corrosion Protection Training Plant

**Product Code :**  
LBNY-0005-11600023



**Description :**

Safety glasses and gloves  
Digital voltmeter  
PC interface for the measurement and data storage  
Digital ammeter on console  
Digital voltmeter on console  
Zn reference cell  
Copper electrode  
Carbon steel electrode  
Cu/CuSO<sub>4</sub> reference cell  
Ag/AgCl reference cell  
Transparent basin to build the electrolytic test bath  
Simple circuit with sliding variable resistor and lamp suitable for inserting into the electrical circuit of the electrolytic cell  
Zinc electrode  
Aluminum electrode  
Magnesium electrode  
DC feeder (with constant voltage, constant current, constant potential facilities)  
Ti-Pt anode  
MMO tubular anode  
Fe bar  
Cu bar  
Fe-Ni bar  
Waterproof resistor with thermostatic device

---

Resistivity fluid cell  
Air pump with relevant sprayer  
Carbon steel electrode (partially coated with epoxy compound)  
Carbon steel electrode (completely coated with epoxy compound)  
Various reagents supplied in plastic cans with safety technical data sheets  
Set of ancillaries and connecting leads  
Set of spare fuses.

**Technical Specification :**

Cathodic Corrosion Protection Training Plant



**Laboratory Instrument India**